

Meetbuisjes

Chromosorb Sorbent Tubes

Sorbent (coating)	Size (mm)	ODxLength	Sections	Sorbent (mg)	Ends	Tube Cover	Separators	Pack Size	Part Number:
Chromosorb 102		8x110	2	100/50	GS	B	WWW	50	226-107
Chromosorb 102		8x110	2	200/100	GS	B	WWW	20	226-106A
Chromosorb 102		8x70	2	100/50	GO	B	WWW	50	226-104 *
Chromosorb 102		6x70	2	66/33	GS	A	WWW	50	226-49-102
Chromosorb 102		6x70	1	100	GS	C	WW	50 Sets	226-105
Chromosorb 106		1/4x3-1/2 in	1	220	GO	B	WW	Each	226-361 ‡ §
Chromosorb 106		1/4x3-1/2 in	1	350	SS	B	SWWS	Each	226-358 ‡ §
Chromosorb 106		10x150	2	600/300	GS	C	WWW	10	226-111A
Chromosorb 106		7x70	2	100/50	GS	B	WWW	50	226-110
Chromosorb 106		6x70	2	75/37	GS	A	WWW	50	226-49-106

Meetapparatuur voor professionals

ProCare | SAFETY
We take intensive Care

www.procaresafety.nl | info@procaresafety.nl

Locatie Groningen

Zernikepark 16a
9747 AN Groningen
T 050 2072172

Locatie Amsterdam

Olympisch Stadion 24-2
1076 DE Amsterdam
T 020 2613999

 @ProCare_safety



Table Notes:

Catalog Number:

- * : Limited shelf-life; contact SKC for more information
- ‡ : Tubes are chemically conditioned before shipping; use within 30 days or recondition. Restocking fee applies.
- § : Each tube has a flow direction arrow and unique number.
- : Limited shelf-life; refrigerator/freezer storage may be required. Contact SKC for more information.

Tube Ends:

- GS: Glass Sealed
- GO: Glass Open
- SS: Stainless Steel Open

Protective Tube Covers:

- A: 6-mm OD x 70-mm length
- B: 8-mm OD x 110-mm length
- C: 10-mm OD x 150-mm length
- D: 10-mm OD x 220-mm length and shorter tubes
- P: Low-volume PUF Tube Holder Do not use an adjustable low flow holder
- V: OVS Tube Holder Do not use an adjustable low flow holder

Separators:

- W: Glass Wool
- G: Glass Fiber Filter
- F: Foam
- N: Nylon Ring
- Q: Quartz Filter
- R: Glass Spacer
- S: Screen
- T: PTFE Ring

Sorbent Equivalencies:

- Anasorb 708 is equivalent to Chromosorb 108.
- Anasorb 727 is equivalent to Chromosorb 106.
- Anasorb 747 is an SKC proprietary sorbent.
- Anasorb C300 is equivalent to Hydrar, Carulite.
- Anasorb CSC is an SKC proprietary sorbent.
- Anasorb GCB1 is equivalent to Carbotrap B (20/40), Carbopack B (60/80).

Trademarks:

- Anasorb ® — SKC Inc
- Chromosorb ® — Manville Corp.
- Porapak ® — Waters Associates, Inc.
- Tenax ® — Buchem B.V.
- XAD ® — Rohm & Haas

Meetapparatuur voor professionals

Locatie Groningen

Zernikepark 16a
9747 AN Groningen
T 050 2072172

Locatie Amsterdam

Olympisch Stadion 24-2
1076 DE Amsterdam
T 020 2613999

